

# USER MANUAL FOR COMPLETE GUIDED SURGERY

---

i bone<sup>®</sup> 

i bone<sup>®</sup> 



**LYRAETK**  
DIGITAL DENTAL PROTOCOLS  
Simple . Safe . Accessible

With research and development experience since 1988 and clinical experience since 1992, **LYRA ETK** has a proven track record in implant design confirmed by the invaluable assistance of internationally recognised research laboratories.

The design of our implants is based on the triple competence of a reactive team experienced in implantology:

- The technical and biomechanical skills of our engineers guarantees the durability of the components and their adaptation to oral applications through modern simulation methods.
- The biological and physiological know-how of the associated laboratories validate the osseointegration capacity of our systems.
- The clinical and practical skills of our dentists and dental technicians ensure the ergonomics of our products, the rationalisation of our protocols and the establishment of ranges adapted to the various clinical cases encountered.

The complete **iBone®** Guided Surgery Kit is also based on the most recent advances in scientific knowledge of implant treatment.

In order to make the most of the complete **iBone®** Guided Surgery Kit, we have produced this manual with a concern for professionalism and we invite you to read it carefully. The smallest detail is important and underscores all the more the difference between the amateur and the specialist.

WARNING	p.4
GENERAL INFORMATION ON COMPLETE GUIDED SURGERY	p.5
Overview of the concept	
Major indications	
Benefits and interests of complete guided surgery	
Limits of complete guided surgery	
PRE-IMPLANTATION STUDY	p.8
Feasibility of implant treatment	
Guide for choosing implants	
DEVELOPMENT OF THE GUIDE	p.11
The GUIDES service on Dental APP	
SURGICAL PREPARATION	p.14
iBone® Complete Guided Surgery Kit	
Positioning the titanium sleeves	
Preparing the surgical kit	
SURGERY - FLAPLESS TECHNIQUE	P.18
Protocol by bone density and implant diameter	
Step-by-step surgical protocol	
Special case: soft tissue that is too thick	
Protecting the connection	
Osseointegration	
In case of failure	

For more information on our products, please visit the websites  
[www.lyraetk.com](http://www.lyraetk.com) et [eshop.lyraetk.com](http://eshop.lyraetk.com)

The e-notice is available online at  
<https://ifu.lyraetk.com>

## WARNING

---

Complete guided surgery is a specific treatment concept, requiring experience with "classical" dental implant dentistry techniques, but also requiring a dual learning:

- Computer-assisted implant planning is a revolutionary tool, but one that requires many clinicians to acquire new computer skills.
- The coordination between you and your team of assistants must be perfect. Indeed, the particular organisation that accompanies this system can be a source of significant preoperative stress for a new user. It is necessary to train your team and give them as many tasks as possible during your interventions so that you can focus on your main objectives.

The experience of the surgeon is crucial, because they must take matters into their own hands in case events do not go as planned.

The instructions and protocols described here must be implemented exclusively using the components and instruments provided by **LYRA ETK**. These instructions will guide you through the different phases of your implant treatments. They are accompanied by the most accurate advice possible, but cannot serve as "recipes", since each clinical case is unique. A very large number of factors act interdependently to achieve a successful implant treatment. It is up to the clinician to understand the key principles, and to draw on their clinical experience.

The technical specifications and clinical advice contained in this manual are for guidance purposes only and shall not give rise to any claim. All the essential information is indicated in the package leaflet provided with the products.

We have taken particular care in the design and manufacture of our products; nevertheless, we reserve the right to make changes or improvements to our complete guided surgery system based on new technical developments. You will be notified of any changes that affect the operating procedure. Depending on the extent of these changes, a new manual may be provided. On the back of your user manual, a version number indicates the date of your publication, which makes it possible to check that you always have the latest updates. You will be able to find the current version of this manual on our website. The reproduction and distribution of all or part of this work requires the prior consent of **LYRA ETK**.

GENERAL INFORMATION  
ON GUIDED SURGERY  
iBONE® COMPLETE GUIDED  
SURGERY

---


### OVERVIEW OF THE CONCEPT

**Complete Guided Surgery (CGS)** can be defined as the ability to transfer a computer-designed implant planning, based on CT scan (CT, Cone Beam) data, into the mouth.

**The transition from a virtual project to a surgical reality is possible** through a surgical guide created by 3D printing from data of the inside of the prosthesis and data of the implant planning made with the planning software

**The advantage of complete guided surgery is undeniable.**

In fact, it offers solutions in the field of pre-implant anatomical investigation (a better predictability of anatomical complications) and thus improves the safety of implant surgery. All important decisions can be taken prior to surgery (choice of implant site, size and diameter of implant, etc.).



### MAJOR INDICATIONS

Complete guided surgery is a treatment applicable to all indications and offers a wide range of flexible solutions:

- From the single case, partial to total edentulism
- Surgical procedures without incision (flapless technique)
- From immediate loading on a frame to deferred loading
- From prefabricated to conventional prosthetic solutions.

**Not to mention that it is indicated when the following criteria are met :**

- The patient meets all the health requirements allowing a surgical intervention to be performed.
- The healing is complete after a graft protocol.
- Satisfactory bone quantity and quality
- Mouth opening sufficient and appropriate **for guided surgical instrumentation.**

## ADVANTAGES AND INTERESTS OF COMPLETE GUIDED SURGERY

- Facilitates decision-making for implant treatment with predictable results, as the system is very reliable, allowing for optimal implant positioning with no margin for error.
- Study of bone density, which will lead to strategic surgical and prosthetic choices: surgical sequence, number of implants, their position, their angulation, type of prosthesis, etc.
- Maximum exploitation of bone volume: no need for timid implantation. Calculation of the volumes to be grafted, which allows the donor site to be chosen.
- Precise orientation of the implant in the mesiodistal and labio-lingual direction. Simplifies protocols and increases the success of surgical treatment.
- Control of bone grafts, or even in some cases, reduction of graft indications by exploiting the remaining bone volume to the maximum.
- Allows perilous or difficult implantations to be performed: search for bicortical supports, lateralisation of the implant path in relation to the dental nerve, etc.
- A further step is taken in the case of flapless surgery, which greatly improves postoperative outcomes, thus saving time and providing substantial comfort for the patient.
- Possibility of combining the osseointegration period and post-extraction implantation in order to have minimum bone resorption with maximum precision.
- Increase patient-practitioner communication.
- This system contributes greatly to ensuring better coordination between the various members of the treatment team (clinician, dental technician, laboratory, assistants), because it has the advantage of being an excellent communication tool and therefore one for teamwork.

## LIMITATIONS OF COMPLETE GUIDED SURGERY

They are mainly technical:

- Insufficient and inadequate mouth opening for **guided surgical instrumentation**.
- The anatomy of the site complicates the insertion of the implant and its anchoring.
- Need for sufficient available bone height.
- The thinness of the bone tables (complete guided surgery requires the existence of an adequate bone quantity and quality).
- Hygiene: the patient must be informed and motivated, but must also have the skills to ensure proper maintenance. Those who are uncooperative will be rejected
- It is important to note that transmucosal surgery is very demanding in terms of indications. There is a high degree of precision with these guides, but it is not absolute and should not overshadow the introduction of sufficient surgical margins.

# PRE-IMPLANTATION STUDY





## PRE-IMPLANTATION STUDY

*It is necessary to assess the feasibility of implant treatment and determine the treatment plan.*

### FEASIBILITY OF IMPLANT

**This study is based on various elements**

- An accurate history based on a medical questionnaire filled in by the patient and collected by the practitioner.
- A clinical and methodical examination of the patient's mouth.
- Biological tests.
- A complete radiological file to determine the available bone volumes.
- Complete study models with both arches in occlusion.
- Implant treatment cannot be initiated until all of the patient's infectious sites have been completely sanitized.

### GUIDE FOR CHOOSING IMPLANTS

**Depending on the amount of bone available**

*In the mesiodistal plane*


- ✓ Allow 2 mm between the coils of an implant and the adjacent natural teeth.
- ✓ Allow 3 mm between the coils of two adjacent implants.

*In the vestibular-palatal-lingual direction*

Leave, if possible, 1.5 to 2 mm of bone thickness around the vestibular, palatal and lingual surfaces.

**According to bone density**


Consider the use of larger implants in low density bones to compensate for the loss of bone/implant contact surface due to cavities.



iBone E & iBone S						
Ø implant body	● Ø 3.5			● Ø 4.3		
Ø threads	3.8	4.3	4.8	4.8	5.5	6.2


# PRE-IMPLANTATION STUDY

## Maxillary




NS	Ø3.5/3.8	RS	Ø4.3/4.8
NM	Ø3.5/4.3	RM	Ø4.3/5.5
NL	Ø3.5/4.8	RL	Ø4.3/6.2

**Legend**



## Mandible



## DEVELOPMENT OF THE GUIDE





## THE « GUIDES » SERVICE

On the Dental APP platform, you can easily order and manage your guides. The guide is developed in several stages.

### ➤ Step 1: Making the impression and the Cone Beam


For the elaboration of the surgical guide, you need two digital images: a CT scan to obtain the anatomical constraints essential to perform the guided surgery (width of the bone crest, positioning of the dental nerve, etc.) and a digital impression to locate the position of future dental implants.

- **You have an intraoral scanner.** Make your impression in the usual way until you obtain the STL file for planning.
- **You do not have an intraoral scanner.** Make your silicone impression in the usual way. The cast model of this impression must then be open-source and CT scanned using a table scanner. Scanning the cast model will provide a STL file for planning.
- Perform the CT scan in a conventional manner as indicated by your equipment. This will provide you with a DICOM file for planning.

### ➤ Step 2 : Implementation of implant planning

- **You have implant planning software.** Import the STL files (impression) and DICOM files (Cone Beam) into the planning software. Perform the planning according to the anatomical constraints of the patient. Once the implant planning has been completed and the surgical guide modelled, you will obtain the STL file necessary for the impression of the guide.

If you want to add the sleeve to your open-source software, the dimension to indicate are follows:



- **You do not have implant planning software.** First, please create your account on the [Dental APP platform](#). If this is already done, then go to your space and click "Create a new case". You must select the level of service that suits you. Lyra Guide can complete the planning and print the guide or one of the two steps according to your choice. Once the service level is selected, follow the instructions on the form. At the end, you will be asked to import the previously obtained STL and DICOM files. Specialists will carry out the implant planning that you will validate before printing.

➤ **Step 3: Printing the surgical guide**

- **You have a 3D printer.** Make an impression of your guide.

*Sleeves are designed to be retained within the guide. However, since each printer is different, we recommend that the printer be set so that the sleeves fit tightly.*

- **You have completed the implant planning but you do not have a 3D printer.** On the Dental APP platform, create a new case and select "Print guide". Follow the form. Import your STL file. LYRA ETK prints your guide.

➤ **Step 4: Placement of the guide sleeves**

- **You have printed the guide yourself.** Insert the sleeves inside the guide. It may be necessary to apply a drop of glue to properly secure the sleeves.
- **LYRA ETK has printed the guide.** The LYRA ETK sleeves are already in place inside the guide. You receive the guide ready to use.

If the case requires the use of anchoring needles, the corresponding sleeves should be placed in the previously defined locations.

*Before starting surgery and after disinfecting the guide, make sure that the surgical guide is properly adapted.*

# SURGICAL PREPARATION



## SURGICAL PREPARATION

### iBONE® COMPLETE GUIDED SURGERY KIT

There are two challenges with making implant sockets:

- ✓ Calibration of the sockets to achieve proper primary stability of the implant, which is essential for osseointegration.
- ✓ The clinician is advised to use sharp drills and to perform the drilling under constant irrigation with refrigerated saline solution. Drills have been developed to allow the saline solution to pass through the guide. It is also important to make back-and-forth movements with the drills to optimise the passage of irrigation within the guide.


The instruments are presented in their order of use from left to right and from top to bottom.



We recommend that you use a "parachute thread" on the instruments to prevent the accidental fall of tools into the patient's throat. It is important to check that the instruments are held correctly on the contra-angle before they are put into the mouth.

Beyond the quality of irrigation, it is also advisable to use drills whose cutting power has not been altered by excessive use.

This kit provides all the necessary tools to complete the guided surgery protocol and to manage all bone densities. It has been designed for the iBone® implant from Ø 3.5/3.8 to Ø 4.3/6.2 and from 6 mm to 14 mm.




1	Adapter		UGC_TE_RP
2	Anchor pins		CGC_AA_15.280
3	Drill for anchor pins		UGC_FA_15.200
4	Gingival punch	TG	TCG_41.49
5		WP	NFE_48.42P
6	Bone drill	TG	FO_47
7		WP	DGCG_41
8	Step drills Ø 2.2	length 6	CGC_NFP_22.060
9		length 8	CGC_NFP_22.080
10		length 10	CGC_NFP_22.100
11		length 12	CGC_NFP_22.120
12	Step drills Ø 2.8	length 14	CGC_NFP_22.140
13		length 6	CGC_NFP_28.060
14		length 8	CGC_NFP_28.080
15		length 10	CGC_NFP_28.100
16	Step drills Ø 3.3	length 12	CGC_NFP_28.120
17		length 14	CGC_NFP_28.140
18		length 6	CGC_NFP_33.060
19		length 8	CGC_NFP_33.080
20	Step drills Ø 3.8	length 10	CGC_NFP_33.100
21		length 12	CGC_NFP_33.120
22		length 14	CGC_NFP_33.140
23	Step drills Ø 4.3	length 6	CGC_NFP_38.060
24		length 8	CGC_NFP_38.080
25		length 10	CGC_NFP_38.100
26		length 12	CGC_NFP_38.120
27	Step drills Ø 4.8	length 14	CGC_NFP_38.140

28	Step drills Ø 4.3	length 6		CGC_NFP_43.060
29		length 8		CGC_NFP_43.080
30		length 10		CGC_NFP_43.100
31		length 12		CGC_NFP_43.120
32		length 14		CGC_NFP_43.140
33	Step drills Ø 4.8	length 6		NFD_38.185
34		length 8		NFD_38.225
35		length 10		NFD_40.185
36	Step drills Ø 5.5	length 12		NFD_40.225
37		length 6		NFV_42.185
38		length 8		NFV_42.225
39		length 10		NFD_45.185
40	Conical drill 3.5 mm	length 12		NFD_45.225
41		length 14		UFE_41.35
42	Conical drill 4.3 mm			UFE_49.42
43	Implant holder tool			CPI_I70
44	Ratchet wrench			CCC_120
45	External hex screwdriver			CCL_HE12.22
46	External hex mandrel			CMA_HE12.22
47	Direct implant holder	TG		CGC_CCP_35
48	Direct implant mandrel			CGC_CMP_35
49	Direct implant holder	WP		CCP_CG_35AWP
50	Direct implant mandrel			CMP_CG_35A.WP
51	Implant holder	TG		CGC_PI_3N
52	Implant holder screw			NPS_VG16.156
53	Implant holder	WP		UGC_PI_UHD48L



## SURGICAL PREPARATION




The drills have been designed to ensure guidance from the start of drilling, regardless of the length of the implant selected. The drills are guided by two elements:

- A cylindrical 10 mm portion provides guidance in the sleeve.
- An active stepped portion whose first part corresponds to the diameter of the previous drill.

### Positioning of THE TITANIUM SLEEVES

The positioning of titanium sleeves on the surgical guide is fixed. Two sleeve references are available :

- TG sleeve : ref.CGC\_TG; the top of the TP sleeve is located 10 mm from the planned implant neck.
- WP sleeve: ref.UGC\_TE\_WP; the top of the WP sleeve is located 9mm from the planned implant neck.



### Using the WP ADAPTER


When inserting an implant with a WP sleeve, the WP adapter must be used to pass drills n°3 to 7 and n°10.

The distance between the surface of the adapter and the neck of the planned implant is 10mm. This means 9mm for the surface of the sleeve + 1mm for the thickness of the adapter.



### Preparation of THE SURGICAL KIT

The surgical kit has 6 locations to position the finned tubes containing the implants. This allows the implants to be positioned in the order desired for the implant insertion to be carried out.






# FLAPLESS SURGICAL TECHNIQUE






## PROTOCOL FOR COMPLETE GUIDED SURGERY by bone density and implant diameter



i bone® - REF. UGCT\_00



## FLAPLESS SURGICAL TECHNIQUE






## FLAPLESS SURGICAL TECHNIQUE





## FLAPLESS SURGICAL TECHNIQUE



# FLAPLESS SURGICAL TECHNIQUE


## STEP-BY-STEP SURGICAL PROTOCOL

### 1. DISINFECTION OF THE SURGICAL GUIDE

The surgical guide does not support autoclave sterilisation, only disinfection is possible. Soak it in chlorhexidine for up to 20 minutes and put it in the mouth quickly.

### 2. CHECK THAT THE SURGICAL GUIDE HAS GOOD STABILITY

Upon receipt, the surgical guide should be tested in the mouth to ensure its suitability to the patient's anatomy. In the case of a dental support guide, you must check with the display windows that the guide is well-supported on the cusps of the adjacent teeth.



### 3. APPLICATION OF LOCAL ANAESTHESIA

Injections of anesthetic products must be made at a distance from the support area of the guide. In fact, the anesthesia stage is a source of mispositioning in relation to the initial project, the local injection of the anesthetic product creates mucosal tumefactions that are sometimes less noticeable, but that can hinder the insertion of the guide with a dental or mucosal support, or promote its placement in an incorrect position, which can lead to significant instability of the guide, in particular the guide with a mucosal support. Analgesia techniques must be mastered to ensure correct and painless guided surgery.



Photo credit: Dr. Ella Bruno (33)

### 4 - PLACEMENT AND STABILISATION OF THE SURGICAL GUIDE

To ensure the success of guided surgery, the surgical guide needs to be stabilized. This can be done in two ways. The method of stabilizing the guide is to be defined during the design.

#### By dental or mucosal support:

The guide is placed and stabilized by the teeth or mucous membrane. Le guide est posé et stabilisé par les dents ou la muqueuse.



#### By anchoring needles :

The stability of the surgical guide can be achieved using anchoring needles. The position of the needles is decided during the implant planning. The locations of the needles are then printed in the surgical guide in order to guide the Ø 1.5 drill for producing the sockets to receive the anchoring needles.



Crédit photo : Dr Ella Bruno (33)





## FLAPLESS SURGICAL TECHNIQUE

### 5 - REMOVAL OF SOFT TISSUES

Set motor speed to 400 rpm and start irrigation. The cutting of soft tissue at the site that will receive the implant is done using a self-guided gingival cutter.

Insert the cutter into the ring before putting the contra-angle into rotation.

Remove tissues with an appropriate instrument (e.g. a curette) and use a bone drill directly.



### 6 – BONE CREST SHAPING

Set motor speed to 400 rpm and start irrigation. Insert the self-guided drill into the sleeve before putting the contra-angle into rotation.

Shape the bone crest to facilitate the guidance of the first drills



### Using the WP ADAPTER

When inserting an implant with a WP sleeve, you need to use the WP adapter to pass the drills n°3 to n°7 and n°10.

The distance between the adapter surface and the planned implant neck is 10mm.




## FLAPLESS SURGICAL TECHNIQUE

### 7 - FIRST DRILL Ø 2.2

The kit contains 3 levels of initial drills of 2.2 mm diameter to ensure guidance of the drills at the start of drilling.


Make sure to follow the arrows indicated on the kit.

For example, for an implant 14 mm long, you need the 2.2 mm diameter, 6 mm long drill, then the 10 mm long drill, followed by the one that is 14 mm long before continuing to drill directly with the 14 mm long stepped drill. Set the motor speed between 1000 and 1200 rpm. Insert the drill into the ring putting the contra-angle into rotation.



Drill to the fixed drill stop with external irrigation. The drill must progress without being forced. If not, it indicates that bone debris cannot be extracted by climbing up the helix. A simple, well-controlled back-and-forth rhythm will allow a more fluid progression of the drill. This action does not require a reversal of motor direction if it is performed at the right time. If the drill is blocked, it can be released in "reverse" mode.

It is also important to make back-and-forth movements to cool the tips of the drills inside the guide. When drilling, check that the bone is bleeding. If not, scrape the bone a bit with a probe or curette to make it bleed so that it can pass through the guide ring without compromising the stability of the surgical guide. If there is no vascularization, it is preferable to close and wait for vascularization.



*The end of the implant does not fit entirely into the tip of the drill hole. The drill hole made is therefore slightly deeper than the length of the implant (1 mm). This helps to avoid any apical compression.*

## FLAPLESS SURGICAL TECHNIQUE

### 8 - ADDITIONAL DRILLING

Use the diagrams on pages 22 to 29 to determine the range of drills corresponding to the diameter of the selected implant and adapt the implant socket to the bone quality of the site.

When drilling, make sure the bone bleeds.

Otherwise, scratch the bone a little to make it bleed. If there is no vascularization, it is preferable to close and wait for vascularization.

Continue the preparation of the site using self-guided drills at the length corresponding to the implant

#### **Partial drilling**

*In densities D2/D3 and D4, the number of the last stepped drill to be used is marked with an asterisk\*. This indicates that the drilling must be partial, i.e. 2 mm shorter than the length of the implant inserted.*

*To do this, use the shorter length drill. For example, for a 10 mm implant, the 8 mm drill will be used for a partial drilling.*

*Attention: for implants of 6 mm in length, the drill indicated with an asterisk\* (indicating partial drilling) is not to be used*

Note 2: the markers allow the embedding of the tap to be viewed with respect to the adjacent anatomical elements and allows the depth of tapping for long implants to be identified.

Note 3: the depth of tapping can be adjusted according to the quality of the bone and the length of the implant.


## FLAPLESS SURGICAL TECHNIQUE

### 9 - IMPLANT INSERTION

The implant can be inserted in three ways:

- Manually with the direct drive wrench
- At the contra-angle with the direct drive wrench
- Manually with the implant holder

This operation must be performed with the utmost care to eliminate the risk of the implant falling out and to ensure that it does not come into contact with any nonsterile components before insertion into the bone site



### 9.A – WITH THE DIRECT DRIVE MANDREL

After opening the sterile tube, connect the end of the wrench or mandrel directly to the implant without removing it from the tube.

Step 1 - Align the hex of the chuck or key with the internal hex of the implant.

Step 2 - To grasp the implant, rotate the chuck or key in the implant connection clockwise until you feel the implant stop rotating in its insert (a device in the insert limits the rotation of the implant during handling).

Step 3 - Insert the mandrel into the implant, applying a slight force on the mandrel to make it retentive on the implant (5N= 500g).

- a. *The feeling of clipping guarantees retention of the implant on the direct socket.*
- b. *The positioning mark is no longer visible, the instrument is correctly inserted in the implant.*
- c. *The positioning mark is visible, the instrument is not aligned or inserted. In this case, return to step 1.*
- d. *The positioning mark is visible, the instrument is not oriented and inserted. In this case, return to step 2.*

Step 4 - With the mandrel properly inserted into the implant, apply a slight counter-clockwise rotation and gently remove the implant from its packaging.

Step 5 - Carry the implant to the recipient site and present it at the well entrance.

*Note: Secure your handling against the risk of the implant falling to the ground or into your mouth.*





For contra-angle placement, we recommend a speed of 15 to 25 rpm to control the implant descent. The contra-angle placement enables the insertion torque of the implant to be measured, and its primary stability to be determined. We recommend inserting the implant at a minimum of 30 N.cm for delayed loading and more than 40 N.cm for immediate or early functional loading. Do not exceed a torque greater than 70 N.cm.



## FLAPLESS SURGICAL TECHNIQUE

### In a D1-D2 bone

For a D1 - D2 bone, it is recommended to finalize the screwing of the implant with a torque spanner in order to ensure the correct insertion of the implant and to check the tightening torque.



### 9.B - WITH THE IMPLANT HOLDER

*Implant gripping in the tube should be done as follows:*

*Step 1 - Place the implant holder in the implant connection.  
The implant holder has hexagonal indexing.*


*Step 2 - Screw the screw inside the implant holder to the implant (5 N.cm) to ensure that the implant holder is held on the implant.*

*Step 3 - Once the implant holder is secured to the implant, place the wrench provided on the implant holder. The wrench clips on the upper part of the implant holder.*

*Step 4 - Move the implant to the recipient site and introduce the implant at the entrance of the drilling socket and then insert the implant using the wrench.*

*We recommend inserting the implant at a minimum of 30 N.cm for delayed loading and more than 40 N.cm for immediate or early functional loading. Do not exceed a torque greater than 70 N.cm.*

*Once the implants are inserted, remove the implant holders by unscrewing their screws beforehand. Implant holders are used to reinforce the proper support of the guide.*




### 10 – REMOVING THE GUIDE

At the end of the procedure, remove the anchoring needles if there are any. Then remove the guide. Excess soft tissue that may interfere with the insertion of the prosthesis can be removed.



### Special case: soft tissue that is too thick

The guide ring required for the complete guided surgery is placed 10 mm from the implant neck. In some cases, the gum thickness is too great and does not allow the ring to be positioned at 10 mm. In this case, it can be placed during implant planning at 12 or 14 mm from the implant neck. It will therefore be necessary to use drills 2 or 4 mm longer. Note: In this case, it will not be possible to use implant holders. Also, for a ring positioned 12 mm from the implant neck, it will not be possible to insert a 14 mm implant, as the drills to be used would then have to be 16 mm drills, which are not available in the complete guided surgery kit. Similarly, for a ring positioned at 14 mm, it will not be possible to insert a 12 mm or 14 mm implant.




**! NB :** When planning, WP sleeves are positioned 9 mm (or 11 or 13 mm) from the implant platform, as the adapter adds 1 mm. Drills n°3 to n°7 and n°10 have an implant length + 10 mm and drills n°8 and n°9, which are passed without an adapter, have an implant length + 9 mm.

### Special case: if you don't want to use the adapter.

You can make 2 guides: one with a TG sleeve and a second one with a WP sleeve

- Use the TG guide for all drills: follow the recommended drilling protocol.
- Place the WP guide
- Then pass drills n°1 and 2 if you need to pass drills n°8 and 9.

Note: Implants with diameters  $\geq 4.8\text{mm}$  should be inserted with a guide fitted with a WP sleeve.



### Protecting the connection

This is ensured:

- Either by iphysio® Profile Designer depending on the tooth and the shape of the prosthetic cradle you wish to obtain. For more information on the iphysio® Profile Designer, go to [iphysio.dental](http://iphysio.dental)
- Or by a healing abutment if you want to work in 1 surgical step. Select the part according to the shape of the prosthetic cradle you wish to obtain and the prosthetic abutment you wish to use afterwards. Manually screw in the chosen abutment using the external hex key and the torque wrench to 10 N.cm.
- Or by a cover screw in the case of a 2-step surgery. The screw is supplied with the implant; it is housed in the cap and will be placed on the implant with the external hex key at 10 N.cm. In this case, the site is sutured, taking care not to pull too much on the soft tissue to avoid an operation.

### Osseointegration

The usual period for obtaining proper osseointegration is:

- 3 months in the mandible
- 6 months in the maxilla due to different bone quality

The clinician must define this period by taking into account bone quality, primary stability of the implant and the prosthetic plan. In some cases, the clinician may decide to connect the prosthetic parts without waiting for osseointegration. However, the clinician must be able to analyse whether the conditions of the clinical case permit immediate loading.

Studies and scientific data show that immediate functional loading on the mandible is possible when the prosthesis is built on 4 implants connected together. Immediate loading is not recommended on a single-unit implant.



## IN CASE OF FAILURE

To place an implant, try to unscrew it with the direct drive wrench. If this solution is not sufficient, please refer to the instructions for the extraction kit. The site may possibly be reimplanted\*, if the patient is fit to receive a new implant, with an implant of a larger diameter, in the event that the insertion of this implant takes place at the same time. To reimplant the site with a smaller diameter implant, it is advisable to wait for the complete healing of the socket.\*\*

*\*It is important to analyse the causes of failure before considering possible re-implantation.*

*\*\*The clinician determines whether it is appropriate to use a filling material.*



EXTRACTION KIT

KDR\_3N

## STUDIES & PUBLICATIONS



**LYRA-ETK**  
DIGITAL DENTAL PROTOCOLS

Simple . Safe . Accessible

LYRA ETK - 726 rue du Général de Gaulle - 74700 Sallanches - France

Service clients France : T. +33 (0)1 76 42 00 60 | International Customer Service: T. +33 (0)1 76 36 14 04  
sales@lyraetk.com | www.lyraetk.com

Medical devices and procedure pack (combination of medical devices) of Class I, in compliance with European Regulation MDR (EU) 2017/745 (CE marking) and IIa and IIb medical devices in compliance with European Regulation MDR (EU) 2017/745 (CE 0459 marking) and/or European Directive 93/42/EEC (CE 0425 marking). Devices not reimbursed by National Insurance; they must be used by oral health professionals. Please read the instructions and user manual or operating technic carefully.

