




















PROSTHETIC GUIDE

Aesthetica+²



CEMENTED PROSTHESIS ON STRAIGHT AND ANGULATED ABUTMENTS

			NP - Ø 4.2	RP - Ø 4.8	WP - Ø 6.5
TEMPORIZATION					
Healing abutments 10 N.cm	h 1.5		ACI_42.47.15	ACI_48.55.15	-
	h 2		-	-	ACI_65.72.20
	h 3		ACI_42.47.30	ACI_48.55.30	ACI_65.72.30
	h 4.5		ACI_42.47.45	ACI_48.55.45	ACI_65.72.45
Cover screws 10 N.cm			ACI_48.35.00		ACI_65.43.00
IMPRESSION					
Pop-in impression copings + screw: APS VG 20 140 5 N.cm			API_42.85	API_48.85	API_65.85
Pick-up impression copings + screw: APV VT 20 154 5 N.cm			APE_T42	APE_T48	APE_T65
Pop-up impression copings 5 N.cm supplied with an impression cap for Pick-up technique: Ø 4.2 : APS TCP 48 55 ● - Ø 4.8 : APS TCP 48 40 ● - Ø 6.5 : APS TCP 65 40 ● + screw: APS VTB 20 174			APU_T42	APU_T48	APU_T65
Implant analogs			ALA_H42	ALA_H48	ALA_H65
Laboratory guide screw			APS_VG20.140		
Scanbody direct on implant + screw: APS VG 20 140 5 N.cm			ETK_AE.NPSB	ETK_AE.RPSB	ETK_AE.WPSB
Scanbody on Esthetibase			ETK_ESTH.LSB		
Bite registration abutment			ETK_AIEP.PO		

CEMENTED PROSTHESIS ON STRAIGHT AND ANGULATED ABUTMENTS






		NP - Ø 4.2		RP - Ø 4.8	WP - Ø 6.5
FINAL RESTORATION					
			APS_PD48.55		APS_PD65.55
			APS_PA48.15		APS_PA65.15
			APS_PA48.20		APS_PA65.20
			APS_PR42.100	APS_PR48.150	APS_PR65.120
			APS_PS48.100		APS_PS65.100
			APS_PCC48.100		APS_PCC65.100
Burn-out sleeves for straight abutments	Rotational		APS_CCC42.55	APS_CCC48.55	APS_CCC65.55
	Non-rotational		APS_CCO42.55	APS_CCO48.55	APS_CCO65.55
Burn-out sleeves for angulated abutments 15° and 20°			APS_CA42.47	APS_CA48.50	APS_CA65.70
Esthetibase interfaces + screw: APS VF20 93 35 N.cm	Rotational		APC_PTC42	APC_PTC48	APC_PTC65
	Non-rotational		APC_PTO42	APC_PTO48	APC_PTO65

CEMENTED PROSTHESIS ON DIRECT CLIP ABUTMENTS

				RP - Ø 4.8			WP - Ø 6.5	
				H 4	H 5.5	H 7	H 4	H 5.5
FINAL ABUTMENT SEATING								
Direct clip abutments 35 N.cm				APS_PP48.40	APS_PP48.55	APS_PP48.70	APS_PP65.40	APS_PP65.55
TEMPORIZATION								
Protection caps	h 5.8		APS_CP48.40	-	-	APS_CP65.40	-	
	h 7.2		-	APS_CP48.55	-	-	APS_CP65.55	
	h 8.8		-	-	APS_CP48.70	-	-	
IMPRESSION								
Snap-on impression copings				APS_TCP48.40	APS_TCP48.55	APS_TCP48.70	APS_TCP65.40	APS_TCP65.55
Snap-on open impression copings for adjusted abutments				APS_TCP48			APS_TCP65	
Direct clip abutment analogs				APS_H48.40	APS_H48.55	APS_H48.70	APS_H65.40	APS_H65.55
FINAL RESTORATION								
Burn-out sleeves	Rotational	h 7		-			APS_BCC65.70	
		h 10		APS_BCC48.100			-	
	Non-rotational	h 7		-			APS_BCO65.70	
		h 10		APS_BCO48.100			-	
DIRECT CLIP KITS								
Single-unit prosthesis				NPS_KIT_NO4	NPS_KIT_NO5	NPS_KIT_NO7	NPS_KIT_WO4	NPS_KIT_WO5
Multi-unit prosthesis				NPS_KIT_NC4	NPS_KIT_NC5	NPS_KIT_NC7	NPS_KIT_WC4	NPS_KIT_WC5


h = supra-implant height

SCREWED PROSTHESIS ON CONOCTA CONICAL ABUTMENTS

			NP - Ø 4.2	RP - Ø 4.8	WP - Ø 6.5
FINAL ABUTMENT SEATING					
Conocta abutments 35 N.cm	h 1.5		APV_P48.15		-
	h 1.8		-		APV_P65.15
TEMPORIZATION					
Conocta protection caps + screw: APV VF 20 56 10 N.cm			APV_CP42	APV_CP48	APV_CP65
Conocta temporary abutments + screw: APV VF 20 56 25 N.cm			APV_PPT42.100	APV_PPT48.100	APV_PPT65.100
IMPRESSION					
Conocta Pick-up impression copings + screw: APE VTC 20 109 5 N.cm			APE_TC42	APE_TC48	APE_TC65
Conocta abutment analogs			ALA_HC42	ALA_HC48	ALA_HC65
Laboratory guide screw			APV_VG20.150		
Scanbody on Conocta abutments + screw: APV VF 20 56 5 N.cm			ETK_AE.NPP	ETK_AE.RPP	ETK_AE.WPP
FINAL RESTORATION					
Conocta burn-out sleeves + screw: APV VF 20 56 25 N.cm	Rotational		APV_BCC42.100	APV_BCC48.100	APV_BCC65.100
	Non-rotational		APV_BCO42.100	APV_BCO48.100	APV_BCO65.100

h = supra-implant height

OVERDENTURE ON O-RING ABUTMENTS

FINAL ABUTMENTS SEATING			
O-Ring abutments 35 N.cm	h 2		APA_OR48.40
	h 4		APA_OR48.60
IMPRESSION			
Pop-in impression taking on O-Ring abutments			
O-Ring abutment analog			OPS_HOBI
FINAL RESTORATION			
O-Ring Supplied with an O-Ring seal of 60 shores			UPA_FOR52
O-Ring seals	50 shores		UPA_JOR50
	60 shores		UPA_JOR60
	70 shores		UPA_JOR70

h = supra-implant height