


PROSTHETIC GUIDE TISSUE LEVEL

Aesthetica+² i bone® 



HEALING SOLUTIONS - TISSUE LEVEL

FOR TISSUE LEVEL IMPLANTS

IPHYSIO® - Delivered with a fixing screw

	 Medical grade titanium.	 15 N.cm	Platform	
			RP	WP
			Supra-implant height (3mm)	ACI_IR ACI_IW

** equivalent to the height of the healing abutment supplied with the implant*

HEALING ABUTMENTS AND COVER SCREWS




HEALING ABUTMENTS

	 The 3 mm high healing abutments are supplied with the implant.	 Titanium	 10 N.cm	Prosthetic profiles	Supra-implant height	References
					1.5	ACI_42.47.15
					3	ACI_42.47.30
					4.5	ACI_42.47.45
					1.5	ACI_48.55.15
					3	ACI_48.55.30
					4.5	ACI_48.55.45
					2	ACI_65.72.20
					3	ACI_65.72.30
					4.5	ACI_65.72.45

COVER SCREW









	 Titanium	Implant diameter	References
		4.2	ACI_48.35.00
		4.8	
		6.5	ACI_65.43.00

CEMENTED PROSTHESIS ON STRAIGHT AND ANGULATED ABUTMENTS




		NP - Ø 4.2	RP - Ø 4.8	WP - Ø 6.5
IMPRESSION				
Pick-up impression copings + screw: APV VT 20 154 5 N.cm		APE_T42	APE_T48	APE_T65
Pop-up impression copings 5 N.cm supplied with an impression cap for Pick-up technique: Ø 4.2 : APS TCP 48 55 ● - Ø 4.8 : APS TCP 48 40 ● - Ø 6.5 : APS TCP 65 40 ● + screw: APS VTB 20 174		APU_T42	APU_T48	APU_T65
Implant analogs		ALA_H42	ALA_H48	ALA_H65
Laboratory guide screw		APS_VG20.140		
Scanbody direct on implant + screw: APS VG 20 140 5 N.cm		ETK_AE.NPSB	ETK_AE.RPSB	ETK_AE.WPSB

h = supra-implant height

NEXT - CEMENTED PROSTHESIS ON STRAIGHT AND ANGULATED ABUTMENTS




		NP - Ø 4.2	RP - Ø 4.8	WP - Ø 6.5	
FINAL RESTORATION					
Straight abutments + screw: APS_VF20.112 35 N.cm			APS_PD48.55		APS_PD65.55
15° angulated abutments + screw : APS_VF20.70 35 N.cm			APS_PA48.15		APS_PA65.15
20° angulated abutments + screw: APS_VF20.70 35 N.cm			APS_PA48.20		APS_PA65.20
Castable abutment Cobalt-chrome base + screw: APS_VF20.70 35 N.cm			APS_PCC48.100		APS_PCC65.100
Burn-out sleeves for straight abutments	Rotational		APS_CCC42.55	APS_CCC48.55	APS_CCC65.55
	Non-rotational		APS_CCO42.55	APS_CCO48.55	APS_CCO65.55
Burn-out sleeves for angulated abutments 15° and 20°			APS_CA42.47	APS_CA48.50	APS_CA65.70
Esthetibase interfaces + screw: APS_ VF20.93 35 N.cm	Rotational		APC_PTC42	APC_PTC48	APC_PTC65
	Non-rotational		APC_PTO42	APC_PTO48	APC_PTO65
Esthetibase interfaces for iPhysio® protocol + screw : APS_ VF20.93 35 N.cm	Rotational		-	APCC_PTC48 (h = 4.8)	APCC_PTC65 (h = 5.1)
	Non-rotational		-	APCC_PTO48 (h = 4.8)	APCC_PTO65 (h = 5.1)

CEMENTED PROSTHESIS ON DIRECT CLIP ABUTMENTS

			RP - Ø 4.8			WP - Ø 6.5	
			Ht 4	Ht 5.5	Ht 7	Ht 4	Ht 5.5
FINAL ABUTMENT SEATING							
Direct clip abutments 35 N.cm			APS_PP48.40	APS_PP48.55	APS_PP48.70	APS_PP65.40	APS_PP65.55
TEMPORIZATION							
Capuchons de protection	h 5.8		APS_CP48.40	-	-	APS_CP65.40	-
	h 7.2		-	APS_CP48.55	-	-	APS_CP65.55
	h 8.8		-	-	APS_CP48.70	-	-
IMPRESSION							
Snap-on impression copings			APS_TCP48.40	APS_TCP48.55	APS_TCP48.70	APS_TCP65.40	APS_TCP65.55
Snap-on open impression copings for adjusted abutments			APS_TCP48			APS_TCP65	
Direct clip abutment analogs			APS_H48.40	APS_H48.55	APS_H48.70	APS_H65.40	APS_H65.55
FINAL RESTORATION							
Burn-out sleeves	Rotational	h 7		-		APS_BCC65.70	
		h 10		APS_BCC48.100		-	
	Non-rotational	h 7		-		APS_BCO65.70	
		h 10		APS_BCO48.100		-	

h = supra-implant height

SCREWED PROSTHESIS ON CONOCTA CONICAL ABUTMENTS

			RP - Ø 4.8	WP - Ø 6.5
FINAL ABUTMENT SEATING				
Conocta abutments 35 N.cm	h 1.5		APV_PI48.15	-
	h 1.8		-	APV_PI65.15
TEMPORIZATION				
Conocta protection caps + screw: APV VF 20 56 10 N.cm			APV_CP48	APV_CP65
Conocta temporary abutments + screw: APV VF 20 56 25 N.cm			APV_PPT48.100	APV_PPT65.100
IMPRESSION				
Conocta Pick-up impression copings + screw: APE VTC 20 109 5 N.cm			APE_TC48	APE_TC65
Conocta abutment analogs			ALA_HC48	ALA_HC65
Laboratory guide screw			APV_VG20.150	
Scanbody on Conocta abutments + screw: APV VF 20 56 5 N.cm			ETK_AE.RPP	ETK_AE.WPP
FINAL RESTORATION				
Conocta burn-out sleeves + screw: APV VF 20 56 25 N.cm	Rotationnelles		APV_BCC48.100	APV_BCC65.100
	Anti-rotationnelles		APV_BCO48.100	APV_BCO65.100

h = supra-implant height

REMOVABLE PROSTHESIS ON O-RING ABUTMENTS

FINAL ABUTMENT SEATING			
O'Ring Abutment 30 N.cm			NPA_OR31.29MALE
EMPREINTE			
Impression on O'Ring Abutment			
O'Ring Analog			OPS_HOBI
RESTAURATION FINALE			
O'Ring Delivered with 60 shores seal			UPA_FOR52
O'Ring seals	50 shores		UPA_JOR50
	60 shores		UPA_JOR60
	70 shores		UPA_JOR70
O'Ring internal hexagonal mandrel			CMO_HI.2526



LYRA Etk - 726 rue du Général de Gaulle - 74700 Sallanches - France

T : +33 (0)1 76 36 14 04 | sales@lyra.dental | www.lyra.dental



Class I, IIa and IIb medical devices (European Directive 93/42/EEC) conforming to current standards and bearing the CE mark (Class I) and CE0459 (Class IIa and IIb). Not reimbursed by the Social Security, they must be used by oral health professionals. Please carefully read the instructions and the user manual. Euroteknika - 726 rue du Général de Gaulle, 74700 Sallanches (France)

PP_TissueLEVEL_ETK_FR_i06_2022-03